
Preliminary Release Information – BETA status



The Business of Science®

Fusion Remote Control API
Using REST
(REpresentational State Transfer)
© Andor Technology 2019



First get Python - https://www.python.org/ftp/python/3.7.2/python-3.7.2-amd64.exe



Page 2

The Business of Science®

- Download installer from Link above
- 2. Run installer in Express mode
- 3. Start Windows Command Prompt
- 4. Run pip install requests as shown
- Note update in yellow below
- 6. Run update of pip

```
Command Prompt
                                                                                                               Microsoft Windows [Version 10.0.14393]
   2016 Microsoft Corporation. All rights reserved.
 \Users\m.browne>pip install requests
 Downloading https://files.pythonhosted.org/packages/7d/e3/20f3d364d6c8e5d2353c72a67778eb189176f08e873c9900e10c0287b84t
 requests-2.21.0-py2.py3-none-any.whl (57kB)
   100%
                                           61kB 707kB/s
 ollecting chardet<3.1.0,>=3.0.2 (from requests)
 Downloading https://files.pythonhosted.org/packages/bc/a9/01ffebfb562e4274b6487b4bb1ddec7ca55ec7510b22e4c51f14098443b8
 chardet-3.0.4-py2.py3-none-any.whl (133kB)
   100%
                                           143kB 1.3MB/s
 ollecting urllib3<1.25,>=1.21.1 (from requests)
 Downloading https://files.pythonhosted.org/packages/62/00/ee1d7de624db8ba7090d1226aebefab96a2c71cd5cfa7629d6ad3f61b79e
 urllib3-1.24.1-py2.py3-none-any.whl (118kB)
                                          122kB 2.3MB/s
 ollecting certifi>=2017.4.17 (from requests)
 Downloading https://files.pythonhosted.org/packages/9f/e0/accfc1b56b57e9750eba272e24c4dddeac86852c2bebd1236674d7887e8a
 certifi-2018.11.29-py2.py3-none-any.whl (154kB)
   100%
                                           163kB 9.2MB/s
collecting idna<2.9,>=2.5 (from requests)
 Downloading https://files.pythonhosted.org/packages/14/2c/cd551d81dbe15200be1cf41cd03869a46fe7226e7450af7a6545bfc474c9
 idna-2.8-py2.py3-none-any.whl (58kB)
   100%
                                           61kB 1.8MB/s
 nstalling collected packages: chardet, urllib3, certifi, idna, requests
uccessfully installed certifi-2018.11.29 chardet-3.0.4 idna-2.8 requests-2.21.0 urllib3-1.24.1
 ou are using pip version 18.1, however version 19.0.3 is available.
 ou should consider upgrading via the 'python -m pip install --upgrade pip' command.
:\Users\m.browne>python -m pip install --upgrade pip
 Downloading https://files.pythonhosted.org/packages/d8/f3/413bab4ff08e1fc4828dfc59996d721917df8e8583ea85385d51125dceff
pip-19.0.3-py2.py3-none-any.whl (1.4MB)
                                           1.4MB 1.2MB/s
 nstalling collected packages: pip
 Found existing installation: pip 18.1
   Uninstalling pip-18.1:
     Successfully uninstalled pip-18.1
 uccessfully installed pip-19.0.3
 :\Users\m.browne>_
```

© Oxford Instruments 2015

If running Windows 7 then the ACL needs to be updated



The Business of Science

ACL – refers to the Access Control List which provides permissions to access files or other resources.

In order to provide permissions to the current user then the following commands must be run from an <u>administrator</u> <u>command prompt.</u> (right-click Command Prompt in the start menu, pick "Run as administrator").

Run the below two commands to allow port 15120, which is the default Fusion port. You'll need to update the "DOMAIN\username" part based on whatever user account will be running Fusion.

netsh http add urlacl url=http://127.0.0.1:15120/ user=DOMAIN\username netsh http add urlacl url=http://localhost:15120/ user=DOMAIN\username


If using just a local user account on the PC, the user account will be "PCNAME\username". If you're not sure, you can check the exact account name from a normal (non-admin) command prompt by:

echo %USERDOMAIN%\%USERNAME%

Switch ON Fusion REST API



- 1. Start Fusion_REST version 2.1.0.34 or later
- 2. Switch to Preferences View
- 3. Select REST API
- 4. Enable REST_API slider to "ON"
- 5. Default Port 15120 should be left at Default
- 6. Click "Back to Imaging"...



Set up Fusion Protocols for REST control




- 1. Ensure correct Channels exist
- 2. Build Protocols
- 3. Check protocols match script names
- 4. Here we have "protocol1" and "protocol2"



Locate Python shortcuts



- 1. Type "Python" into Windows search box
- 2. Windows results opens Python on top
- 3. Right mouse Open file location click
- 4. Windows File Explorer shows folder
- 5. Double Click on "IDLE (Python 3.7 64 bit)
- 6. Right Mouse "Pin to Taskbar" if required




Load Fusion Python Scripts



X

The Business of Science®

- 1. IDLE editor will Open as shown
- 2. Locate Fusion REST Docs and Examples



Python 3.7.2 Shell


Page 7

Edit Fusion Python Scripts



The Business of Science®

- 1. mit-test.py opens in its own window
- 2. Edit protocol names as required
- 3. Change dummy processes to real triggers etc



Page 8

Run Fusion Python Scripts



- 1. Select "Run" from the menu
- 2. Run Module mit-test.py
- 3. Observe Fusion executes scripts, IDLE reports output

```
Python 3.7.2 Shell
                                                                                 X
File Edit Shell Debug Options Window Help
Python 3.7.2 (tags/v3.7.2:9a3ffc0492, Dec 23 2018, 23:09:28) [MSC v.1916 64 bit
(AMD64)] on win32
Type "help", "copyright", "credits" or "license()" for more information.
 RESTART: C:\Users\m.browne\AppData\Local\Programs\Python\Python37\Tools\scripts
\mit-test.pv
initialising
==== start loop #1 ====
(pretending to wait for input trigger)
== run protocol1 ==
(pretending to send output trigger)
== run protocol2 ==
(pretending to send output trigger)
==== start loop #2 ====
(pretending to wait for input trigger)
== run protocol1 ==
(pretending to send output trigger)
== run protocol2 ==
(pretending to send output trigger)
==== start loop #3 ====
(pretending to wait for input trigger)
== run protocol1 ==
(pretending to send output trigger)
== run protocol2 ==
(pretending to send output trigger)
success
>>>
```

End of Preliminary setup



- 1. Now you can create your own protocol sequences
- 2. Control external hardware through Python or other
- 3. Note that fusionrest.py includes only the definitions
- 4. This means they do not execute but define REST methods
- 5. For more help contact the Andor team
- 6. email: TBD
- 7. The Fusion team will be engaged as needed.

```
fusionrest.py - C:\Users\m.browne\AppData\Local\Programs\Python\Python37\Tools\scripts\fusionrest.py (3.7.2)
                                                                                                              File Edit Format Run Options Window Help
import requests
import json
import time
host = "localhost"
port = 15120
class ApiError (Exception):
        Indicates an error while calling the Fusion REST API.
        def __init__(self, endpoint, code, reason):
                Creates an new 'ApiError' instance.
                self. endpoint = endpoint
                self. code = code
                self. reason = reason
                return "<ApiError at {}: {} {}>".format(self. endpoint, self. code, self. reason)
        def __str__(self):
                return self.__repr__()
        def endpoint(self):
                 Gives the name of the API endpoint for which the error happened.
                return self. endpoint
        def code(self):
                Gives the HTTP response code for the error, as returned by the API.
                Also see `.reason()` for a more readable description of the problem.
                return self. code
        def reason(self):
                 Gives the reason for the error, as returned by the API. (a string)
                return self. reason
def make address(endpoint):
        return "http://{}:{}{}" format(host. nort. endnoint)
                                                                                                               Ln: 1 Col: 0
```



The Business of Science®

© Oxford Instruments 2015